# **Yukon<sup>®</sup> Chrome PC**

---

Sirolimus Eluting Cobalt Chromium Coronary Stent System  
Endoprothèse coronaire (stent) en cobalt-chrome à libération de sirolimus  
Sistema di stent coronarico a eluizione di cromo-cobalto Sirolimus  
Sistema de stents coronarios de cromo cobalto liberadores de Sirolimus  
Sistema de Stent Coronário em Cromo-cobalto com Eluição de Sirolimus

---

**INSTRUCTIONS FOR USE**

**MODE D'EMPLOI**

**ISTRUZIONI D'USO**

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## Table of Content

1. Product Description .....	3
2. Indications .....	5
3. Contra-indications .....	5
4. Warning .....	5
5. Precautions .....	5
6. Adverse Events .....	8
7. Patient Counselling Information.....	9
8. Packaging .....	9
9. Operator's Manual .....	9
10. Disclaimer of Warranty and Limitation of Remedy .....	11

These instructions are intended for all personnel who will handle the „Yukon® Chrome PC“, a premounted stent system, and should be read carefully before the product is used!

## Product Label – Symbol Description

	Nominal length of stent		Temperature limit
	Rated burst pressure		Nominal diameter of stent
	Catalogue number		Nominal pressure
	CE Mark		Batch Code
	Manufacturer		MR conditional
	Date of manufacture		Expiry Date
	Consult instructions for use		Caution
	Do not re-use		Do not use if packaging is damaged
	Do not resterilize		Keep away from sunlight
	Keep dry		Sterilized using ethylene oxide

## 1. Product Description

The **Yukon Chrome PC Sirolimus Eluting Coronary Stent system** is a combination product comprised of two regulated components: A Device (A Stent System) and a drug product (A formulation of Sirolimus in Resomer – R202S).

### 1.1 Device Component Description

The device components consist of a stent mounted on a Stent delivery System. The Physical characteristics of the device components are shown below:

Brand Name	Yukon Chrome PC Sirolimus Eluting Coronary Stent System
Available Stent Lengths, Unexpanded (mm):	8, 12, 16, 18, 21, 24, 28, 32, 40
Available Stent Diameters (mm):	2.0, 2.25, 2.5, 2.75, 3.0, 3.5, 4.0
Stent Material:	Cobalt Chromium Alloy, laser cut from seamless tubing in a serpentine pattern coated with a Biodegradable Polymer and Sirolimus Combination.
Strut Thickness:	68 µm / 79 µm / SV / MV Design
Drug	Sirolimus
Nominal Stent Foreshortening:	None
Delivery System Usable Length:	143 cm
Delivery System Y-Adapter Ports:	Single access port to inflation/deflation lumen. A guide wire exit port is located at 25 cm ... 29 cm from the tip. Designed for guide wire 0.014" (0.36 mm)
Stent Delivery Balloon	Polyamide balloon, nominally equal to stent. Mounted stent length and location is defined by radio opaque markers proximal and distal to stent.
Balloon Inflation Pressure	Nominal Inflation Pressure: 1100 kPa (11 bar) Rated Burst Pressure: 1600 kPa (16 bar)
Guiding Catheter Inner Diameter	5 F (0.065" ID) for all diameters
Catheter Shaft Outer Diameter	1.9F proximally, 2.7F distally
Guide Wire	0.014 inch

### 1.2 Drug Component Description


A brief description of the drug and the therapeutic class to which it belongs.

#### 1.2.1 Sirolimus

Pharmacologic class: Macrocyclic lactone

Therapeutic class: Immunosuppressant

Pregnancy risk: Category C

(3S, 6R, 7E, 9R, 10R, 14S, 15E, 17E, 19E, 21S, 23S, 26R, 27R, 34aS)-9, 10, 12, 13, 14, 21, 22, 23, 24, 25, 26, 27, 32, 33, 34, 34a-hexadecahydro-9, 27-dihydroxy-3-[[[(1R)-2-[(1S, 3R, 4R)-4-hydroxy-3-methoxycyclohexyl]-1-methyl]-10, 21-dimethoxy-6, 8, 12, 14, 20, 26-hexamethyl-23, 27-epoxy-3H-pyrido[2, 1-c][1, 4]oxazacyclohentacontine-1, 5, 11, 28, 29(4H, 6H, 31H)-pentone]

International Non-proprietary name: Rapamycin (Sirolimus)

Nationally Approved Name: Sirolimus

Molecular Formula: C<sub>51</sub>H<sub>79</sub>NO<sub>13</sub>

Molecular Weight: 914.19 g/mol

Sirolimus is off-white crystalline powder, highly lipophilic, insoluble in water, sparingly soluble in Methanol, highly soluble in dichloromethane.

### 1.2.1.2 Yukon Chrome PC Nominal Sirolimus dosages as per various lengths

**Yukon Chrome PC** Sirolimus Eluting Coronary Stent has Sirolimus as the active pharmaceutical ingredient which releases from the matrix of Resomer 202S. The active ingredient, Sirolimus loading per stent is given in the table below:

Ref No.	Nominal unexpanded Stent Length (mm)	Nominal Stent Diameter at NP (mm)	Nominal Sirolimus content (µg)	Ref No.	Nominal unexpanded Stent Length (mm)	Nominal Stent Diameter at NP (mm)	Nominal Sirolimus content (µg)
T-CMG2008PC	8	2.0	100	T-CMG2521PC	21	2.5	263
T-CMG22508PC	8	2.25	100	T-CMG27521PC	21	2.75	263
T-CMG2508PC	8	2.5	100	T-CMG3021PC	21	3.0	263
T-CMG27508PC	8	2.75	100	T-CMG3521PC	21	3.5	263
T-CMG3008PC	8	3.0	100	T-CMG4021PC	21	4.0	263
T-CMG3508PC	8	3.5	100	T-CMG2024PC	24	2.0	300
T-CMG4008PC	8	4.0	100	T-CMG22524PC	24	2.25	300
T-CMG2012PC	12	2.0	150	T-CMG2524PC	24	2.5	300
T-CMG22512PC	12	2.25	150	T-CMG27524PC	24	2.75	300
T-CMG2512PC	12	2.5	150	T-CMG3024PC	24	3.0	300
T-CMG27512PC	12	2.75	150	T-CMG3524PC	24	3.5	300
T-CMG3012PC	12	3.0	150	T-CMG4024PC	24	4.0	300
T-CMG3512PC	12	3.5	150	T-CMG2028PC	28	2.0	350
T-CMG4012PC	12	4.0	150	T-CMG22528PC	28	2.25	350
T-CMG2016PC	16	2.0	200	T-CMG2528PC	28	2.5	350
T-CMG22516PC	16	2.25	200	T-CMG27528PC	28	2.75	350
T-CMG2516PC	16	2.5	200	T-CMG3028PC	28	3.0	350
T-CMG27516PC	16	2.75	200	T-CMG3528PC	28	3.5	350
T-CMG3016PC	16	3.0	200	T-CMG4028PC	28	4.0	350
T-CMG3516PC	16	3.5	200	T-CMG2032PC	32	2.0	400
T-CMG4016PC	16	4.0	200	T-CMG22532PC	32	2.25	400
T-CMG2018PC	18	2.0	225	T-CMG2532PC	32	2.5	400
T-CMG22518PC	18	2.25	225	T-CMG27532PC	32	2.75	400
T-CMG2518PC	18	2.5	225	T-CMG3032PC	32	3.0	400
T-CMG27518PC	18	2.75	225	T-CMG3532PC	32	3.5	400
T-CMG3018PC	18	3.0	225	T-CMG4032PC	32	4.0	400
T-CMG3518PC	18	3.5	225	T-CMG27540PC	40	2.75	500
T-CMG4018PC	18	4.0	225	T-CMG3040PC	40	3.0	500
T-CMG2021PC	21	2.0	263	T-CMG3540PC	40	3.5	500
T-CMG22521PC	21	2.25	263	T-CMG4040PC	40	4.0	500

### 1.2.2 Mechanism of action

It is known that sirolimus inhibits T-lymphocyte activation and smooth muscle and endothelial cell proliferation in response to cytokine and growth factor stimulation. In cells, sirolimus binds to the immunophilin, FK Binding Protein-12 (FKBP-12). The sirolimus-FKBP-12 complex binds to and inhibits the activation of the mammalian Target of Rapamycin (mTOR), leading to inhibition of cell cycle progression from the G1 to the S phase.

## 2. Indications

The **Yukon Chrome PC** Sirolimus Eluting Stent is indicated for improving luminal diameter and reducing restenosis for treatment of coronary artery lesions in native coronary arteries ranging from 2.0 mm to 4.0 mm.

The **Yukon Chrome PC** Sirolimus Eluting Cobalt Chromium Coronary Stent System is indicated for intraluminal chronic placement in stenosed coronary artery or aorto-coronary bypass grafts in order to obtain vessel patency following acute or subacute coronary artery obstruction. It is also indicated for restenosis or arterial dissection after PTCA procedure. Patients being considered for stent implantation should be acceptable candidates for coronary balloon angioplasty.

## 3. Contra-indications

Use of the **Yukon Chrome PC** Sirolimus Eluting Coronary Stent System is contraindicated in the following patient types:

- Patients with a contraindication for anti-platelet/anti-coagulant therapy.
- Patients in whom aortocoronary bypass procedures are contraindicated.
- Patients with imminent thrombus formation.
- Patients with proximal atherosclerosis and/or significant vessel tortuosity in whom guide wire access and adequate guiding catheter support is prohibited.
- Coronary artery spasm in the absence of a significant stenosis.
- Patients with unprotected left main coronary vessel disease.
- Target lesion distal to previously placed stents.
- Lesions that cannot be successfully pre-dilated.
- Patients judged to have lesion that prevents complete inflation of an angioplasty balloon.
- Known hypersensitivity to Sirolimus or its derivatives.
- Known hypersensitivity to Poly-lactide PLA (Polymers might enhance inflammatory reactions and prothrombotic response).
- Known allergy to Cobalt Chromium Alloy.
- Patients in whom the necessary pharmaceutical treatments are contraindicated.
- Severe allergic reaction to the contrast medium, stent and/or delivery system material and/or the necessary pharmaceutical treatment (e.g. aspirin).
- Pregnant women and women of childbearing potential.

## 4. Warning

- Please ensure that the inner package has not been opened or damaged as this may indicate the sterile barrier has been breached.
- Since the use of this device carries the associated risk of sub acute thrombosis, vascular complications and/or bleeding events, judicious selection of patients is necessary.
- Persons allergic to Cobalt Chromium Alloy or Sirolimus may suffer an allergic reaction to this implant.
- Do not use Drug eluting Balloon before Drug Eluting Stent implantation.
- In case of serious pre-existing diseases/conditions (e.g. renal failure) or in case of uncooperative patients surgeon should carefully consider and assess additional risks for the patient before using the product.

## 5. Precautions

### 5.1 General Precautions

1. The **Yukon Chrome PC** Sirolimus eluting stent should not be exposed to any direct handling or contact with liquids prior to preparation and delivery as the coating may be susceptible to damage or premature drug elution.
2. Only physicians who have received adequate training should perform implantation of the stent.
3. Stent placement should only be performed at hospitals where emergency coronary artery bypass graft surgery can be readily performed.
4. Subsequent stent blockage may require repeat dilatation of the arterial segment containing the stent. The long-term outcome following repeat dilatation of endothelialised stents is not well characterized.

### 5.2 Use of Multiple Stents

The extent of the patient's exposure to drug and polymer is directly related to the number of stents implanted. Use of more than two **Yukon Chrome PC** Stents has not been fully evaluated. Use of more than two **Yukon Chrome PC** Stents will result in the patient receiving larger amounts of drug. When multiple stents are required, resulting in stent-to-stent contact, stent materials should be same type of material to avoid the possibility of dissimilar metal corrosion.

## 5.3 Brachytherapy

The safety and effectiveness of the **Yukon Chrome PC Stent** in patients with prior brachytherapy of the target lesion have not been established. The safety and effectiveness of use of brachytherapy to treat in-stent restenosis in **Yukon Chrome PC Stent** have not been established.

## 5.4 Use in Conjunction with Other Procedures

The safety and effectiveness of using mechanical atherectomy devices (directional atherectomy catheters, rotational atherectomy catheters) or laser angioplasty catheters in conjunction with **Yukon Chrome PC Stent** implantation have not been established.

## 5.5 Use in Special Populations

### 5.5.1 Pregnancy

There are no adequate and well-controlled studies in pregnant women or men intending to father children. Systemic levels of Sirolimus have not been demonstrated in any pre-clinical or clinical trials with the **Yukon Chrome PC Stent**. Effective contraception should be initiated before implanting an **Yukon Chrome PC Stent** and for 12 weeks after implantation. The **Yukon Chrome PC Stent** should be used during pregnancy only if the potential benefit outweighs the potential risk to the embryo or fetus.

### 5.5.2 Use during Lactation

A decision should be made whether to discontinue nursing or to implant the stent, taking into account the importance of the stent to the mother.

Sirolimus is excreted in trace amounts in milk of lactating rats. It is not known whether sirolimus is excreted in human milk. The pharmacokinetic and safety profiles of sirolimus in infants are not known. Because many drugs are excreted in human milk and because of the potential for adverse reactions in nursing infants from sirolimus, a decision should be made whether to discontinue nursing or to implant the stent, taking into account the importance of the stent to the mother.

### 5.5.3 Pediatric Use

The safety and efficacy of the **Yukon Chrome PC Stent** in pediatric patients have not been established.

### 5.5.4 Geriatric Use

Clinical studies of the **Yukon Chrome PC Stent** did not find that patients age 65 years and over differed with regard to safety and efficacy compared to younger patients.

## 5.6 Lesion/Vessel Characteristics

The safety and effectiveness of the **Yukon Chrome PC Sirolimus Eluting Coronary Stent** have not been established in the following patient populations:

1. Patients with unresolved vessel thrombus at the lesion site.
2. Patients with coronary artery reference vessel diameter < 2.0 mm or >4.0 mm.
3. Patients with lesions located in the left main coronary artery, ostial lesions, or lesions located at a bifurcation.
4. Patients with diffuse disease or poor overflow distal to the identified lesions.
5. Patients with tortuous vessels in the region of the obstruction or proximal to the lesion.
6. Patients with a recent acute myocardial infarction where there is evidence of thrombus or poor flow.

**In general, for implantation of drug-eluting stents in different lesions and patient settings the latest clinical guidelines, e.g. ESC guidelines for myocardial revascularization 2018, should be considered.**

## 5.7 Drug Interactions

Consideration should be given to the potential for drug interaction when deciding to place the **Yukon Chrome PC Stent** in a patient who is taking a drug that could interact with Sirolimus, or when deciding to initiate therapy with such a drug in a patient who had recently received the **Yukon Chrome PC Stent**. The drug interactions are unlikely to be detectable.

## 5.8 Magnetic Resonance Imaging (MRI) – Stent Migration

Non-clinical testing has demonstrated that **Yukon Chrome PC Coronary Stent**, in single configuration up to 40 mm (nominal), is MR Conditional. A patient with a **Yukon Chrome PC Coronary Stent** can be scanned safely, immediately after placement of this implant, under the following conditions:

Static magnetic field of 3 Tesla, 1,5 Tesla or less

Spatial gradient field of 10 T/m or less

Maximum whole-body-averaged specific absorption rate (SAR) of 2,0 W/kg (under normal operating mode only) for 15 minutes of scanning or less. The **Yukon Chrome PC** Coronary Stent should not migrate in this MRI environment.

Stent heating was derived by relating the measured non-clinical, in vitro temperature rise in a 3 Tesla MR-scanner and in a 1,5 Tesla MR-scanner to the local specific absorption rates (SARs). Tests were performed in a 3 Tesla MR-scanner Skyra (Siemens (Erlangen, Germany) Software: Numaris/4D13C) and in a 1,5 Tesla MR-scanner Vision (Siemens (Erlangen, Germany), Software: Numaris/4, A35). The maximum whole body averaged SAR was determined by calculation. At a length of up to 41 mm, a single stent produced a non-clinical maximum local temperature rise of 3,0 K at a maximum whole body averaged SAR of 2,0 W/kg (normal operating mode) for 15 minutes. These calculations do not take into consideration the cooling effects of blood flow.

Non-clinical testing at field strengths greater than 3 Tesla has not been performed to evaluate stent migration or heating. Non-clinical testing was performed with a single stent. Effects produced by overlapping stents or on stents with fractured struts, especially on heating in the MRI environment, have not been tested and can be significantly higher.

As demonstrated in non-clinical testing, an image artifact can be present when scanning the **Yukon Chrome PC** Coronary Stent. MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the stent. Therefore it may be necessary to optimize the MR imaging parameters for the presence of **Yukon Chrome PC** Coronary Stents.

## 5.9 Stent Handling Precautions

1. For single use only. Do not resterilize or reuse this device. Note the "Use Before" date on the product label.
2. Do not remove the stent from the delivery balloon as removal may damage the stent and/or lead to stent embolization. The stent system is intended to perform as a system.
3. Do not induce a vacuum on the delivery system prior to reaching the target lesion.
4. Special care must be taken not to handle or in any way disrupt the stent on the balloon. This is most important during stent system removal from packaging, placement over guide wire, and advancement through rotating haemostatic valve adaptor and guiding catheter hub.
5. Do not "roll" the mounted stent with your fingers as this action may loosen the stent from the delivery balloon and may damage the coating.
6. Use only the appropriate balloon inflation media. Do not use air or any gaseous medium to inflate the balloon as this may cause uneven expansion and difficulty in deployment of the stent.
7. In the event the **Yukon Chrome PC** Stent is not deployed, follow product return procedures and avoid handling of the stent with hands.

## 5.10 Stent Placement Precautions

- Do not prepare or pre-inflate balloon prior to stent deployment other than as directed.
- When treating multiple lesions, the distal lesion should be initially stented, followed by stenting of the proximal lesion. Stenting in this order obviates the need to cross the proximal stent in placement of the distal stent and reduces the chances for dislodging the proximal stent. Implanting a stent may lead to dissection of the vessel distal and/or proximal to the stent and may cause acute closure of the vessel requiring additional intervention (CABG, further dilatation, placement of additional stents, or other). Do not expand the stent if it is not properly positioned in the vessel (see 5.11 Stent/System Removal Precautions). Placement of a stent has the potential to compromise side branch patency. The vessel should be pre-dilated with an appropriate sized balloon.
- Balloon pressures should be monitored during inflation. Do not exceed rated burst pressure as indicated on the product label. Use of pressures higher than those specified on the product label may result in a ruptured balloon with possible intimal damage and dissection.
- Do not attempt to pull an unexpanded stent back through the guiding catheter, as dislodgement of the stent from the balloon may occur. Remove as a single unit per instructions in section 5.11 (Stent/System Removal Precautions).
- An unexpanded stent should be introduced into the coronary arteries one time only. An unexpanded stent should not be subsequently moved in and out through the distal end of the guiding catheter as stent damage or stent dislodgement from the balloon may occur.
- Stent retrieval methods (use of additional wires, snares and/or forceps) may result in additional trauma to the coronary vasculature and/or the vascular access site. Complications may include bleeding, hematoma or Pseudoaneurysm.

## 5.11 Stent/System Removal Precautions

- Should unusual resistance be felt at any time during either lesion access or removal of the Delivery System post-stent implantation, the entire system should be removed as a single unit.
- Do not attempt to pull an unexpanded stent back through the guiding catheter while engaged in the coronary arteries, as stent damage or stent dislodgement from the balloon may occur.

**When removing the Delivery System as a single unit**

- DO NOT retract the Delivery System into the guiding catheter.
- Position the proximal balloon marker just distal to the tip of the guiding catheter.
- Advance the guide wire into the coronary anatomy as far distally as safely possible.
- Tighten the rotating haemostatic valve to secure the Delivery System to the guiding catheter. Then remove the guiding catheter, guiding wire and Delivery System as a single unit.

Failure to follow these steps and/or applying excessive force to the Delivery System can potentially result in loss or damage to the stent and/or Delivery System components.

If it is necessary to retain guide wire position for subsequent artery/lesion access, leave the guide wire in place and remove all other system components.

**5.12 Post Implantation Precautions**

- Great care must be exercised when crossing a newly deployed stent with a coronary guide wire or balloon catheter to avoid disrupting the stent geometry and stent coating.
- Prescribe an antiplatelet therapy for a period of 6 months to reduce the risk of stent thrombosis and follow the latest clinical guidelines for antithrombotic therapy, e.g. ESC guidelines 2018.

**6. Adverse Events****6.1 Potential Adverse Events**

Potential adverse events which may be associated with the use of a coronary stent include but are not limited to:

- Abrupt stent closure
- Acute myocardial infarction
- Allergic reaction to anti-coagulant and/or antithrombotic therapy or contrast medium
- Angina
- Aneurysm
- Arrhythmias, including ventricular fibrillation (VF) and ventricular Tachycardia (VT)
- Arterial perforation
- Arterial rupture
- Arteriovenous fistula
- Bleeding complications
- Bradycardia
- Cardio Tamponade
- Cardiogenic Shock
- Coronary spasm
- Coronary or stent embolism
- Coronary or stent thrombosis
- Death
- Dissection of the coronary artery
- Drug reactions to antiplatelet agents / anticoagulation agents / contrast medium
- Emboli, distal (air, tissue or thrombotic emboli)
- Emergency or non-emergent Coronary Artery Bypass Graft Surgery
- Entry site complications
- Heart Failure
- Hematoma
- Hemorrhage, requiring transfusion
- Hypotension / Hypertension
- Infection
- Infection and/or pain at the access site
- Injury to the coronary artery
- Ischemia
- Nausea and vomiting
- Palpitations
- Perforation or rupture
- Pericardial effusion
- Pseudoaneurysm, femoral
- Renal Failure
- Respiratory Failure
- Restenosis of the stented segment

- Rhythmical disturbances
- Shock/Pulmonary edema
- Spasm
- Stroke/cerebrovascular accident/TIA
- Total occlusion of the coronary artery
- Unstable angina pectoris
- Vascular complications, which may require vessel repair
- Ventricular fibrillation

Potential adverse events not captured above, that may be unique to the Sirolimus drug coating:

- Allergic/immunologic reaction to drug or stent coating
- Alopecia
- Anemia
- Blood product transfusion
- Gastrointestinal symptoms
- Hematologic dyscrasia (including leukopenia, neutropenia, thrombocytopenia)
- Hepatic enzyme changes
- Histologic changes in vessel wall, including inflammation, cellular damage or necrosis
- Myalgia/Arthralgia
- Peripheral neuropathy

## 7. Patient Counselling Information

Physicians should consider the following in counselling patient about this product:

- Discuss the risks associated with stent placement
- Discuss the risks associated with a Sirolimus Eluting implant
- Discuss the risks/benefits issues for this particular patient
- Discuss alteration to current lifestyle immediately following the procedure and over the long terms.

## 8. Packaging

- Content: One (1) Yukon Chrome PC Sirolimus Eluting Cobalt Chromium Stent System.
- Sterile: This device is sterilized with ETO, Non-pyrogenic.
- Do not use if the package is opened or damaged.
- Storage: Store in cool and dry place at temperature between 8°C to 25°C.

## 9. Operator's Manual

### 9.1 Access to Package Holding Sterile Stent Delivery System

- Tear open outer foil pouch to reveal second inner pouch.  
Note: DO NOT drop or hand inner pouch into sterile field.
- Remove inner pouch from outer foil pouch.
- Peel open inner pouch using aseptic technique to reveal sterile product.

### 9.2 Inspection Prior to Use

Do not use if there is any damage to the packaging. Prior to using the **Yukon Chrome PC** Stent System, carefully remove the system from the package and inspect for bends, kinks, and other damage. Verify that the stent is located between the radiopaque balloon markers.

### 9.3 Materials Required

Material
<ul style="list-style-type: none"> <li>• Introducer sheath compatible with the guiding catheter</li> <li>• Femoral or brachial guiding catheter in the appropriate size and configuration to match the coronary artery (between 5F - 9F).</li> <li>• Haemostatic valve</li> <li>• 60% contrast medium diluted 1:1 with heparinized physiological saline</li> <li>• 10 - 20 ml Luer-Lock syringe (optional)</li> <li>• Inflation device with manometer</li> <li>• Guide wire 0.014" (0.36mm)</li> <li>• Guide wire introducer (optional)</li> <li>• Guidewire torque device (optional)</li> </ul>

## 9.4 Preparation

### 9.4.1 Guide wire Lumen Flush

Step	Action
1.	Prepare the inflation device according to the manufacturer's instructions with diluted contrast medium.
2.	Prepare diluted contrast medium and sterile saline solution.
3.	Fill a 10 or 20 cc syringe with 4 cc of contrast medium.
4.	Carefully remove the stent system from the sterile packaging. Use usual sterile technique to do this. Pull the stent system out of the transport tubing. Remove the transport protection. CAUTION: Do not manipulate, touch or handle the stent with your fingers as this may cause contamination or dislodgement of the stent from the delivery balloon.
5.	Remove the product protective stent sheath by grasping the sheath at the distal end and gently remove distally. If unusual resistance is felt during product mandrel and stent sheath removal, do not use the product and replace with another.

### 9.4.2 Delivery System Preparation

Step	Action
1.	Attach the inflation device/syringe to the stopcock. Attach to inflation port.
2.	With tip down, orient Delivery System vertically.
3.	Open stopcock to Delivery System. Pull negative for 30 seconds. Release to neutral for contrast fill.
4.	Close stopcock to Delivery System. Purge inflation device/syringe of all air.
5.	Repeat steps 3 to 5 until all air is expelled. Note: If air is seen in shaft, repeat balloon preparation steps 3 to 5 to prevent uneven stent expansion.
6.	If a syringe was used, attach a prepared inflation device to stopcock.
7.	Open stopcock to Delivery System.
8.	Leave on neutral.

## 9.5 Delivery Procedure

Step	Action
1.	Prepare the vascular access site according to standard practice.
2.	Preciliate the lesion with a PTCA catheter.
3.	Maintain neutral pressure on the inflation device. Open the rotating hemostatic valve as widely as possible.
4.	Backload the delivery system onto the proximal portion of the guide wire while maintaining the guide wire position across the target lesion.
5.	Advance the stent delivery system over the guide wire to the target lesion. Use the radiopaque balloon markers to position the stent across the lesion; perform angiography to confirm the position of the stent. NOTE: If during the process of moving the Delivery System into position you notice the stent has moved on the balloon, do not deploy the stent. The entire system should be removed as a single unit (see section 5.11 (Stent/System Removal Precautions) for specific Delivery System removal instructions).
6.	Tighten rotating hemostatic valve. Stent is now ready to be deployed.

## 9.6 Deployment Procedure

Step	Action
1.	CAUTION: Refer to product label for in vitro stent diameter, deployment pressure, and RBP. Deploy the stent by slowly pressurizing the system in 2 bar increments, every 5 seconds, until the stent is completely expanded. Maintain pressure for 45 seconds after inflation of stent. If necessary, the balloon can be repressurized or further pressurized to assure complete apposition of the stent to artery wall. Do not exceed RBP.
2.	Deflate balloon by pulling negative on inflation device for 30 seconds.
3.	Await for removal until the balloon is completely evacuated.

## 9.7 Removal Procedure

Step	Action
1.	Ensure that the balloon is fully deflated.
2.	Fully open rotating haemostatic valve.
3.	While maintaining guide wire position and negative pressure on inflation device, withdraw Delivery System. NOTE: Should unusual resistance be felt at any time during either lesion access or removal of Delivery System post-stent implantation, the entire system should be removed as a single unit (see section 5.11 (Stent/System Removal Precautions) for specific Delivery System removal instructions).
4.	Tighten rotating haemostatic valve.
5.	Repeat angiography to assess stented area. If necessary, post dilates within stent. Balloon inflations should incorporate balloon size closely matching vessel.
6.	Final stent diameter should match reference vessel. ASSURE STENT IS NOT UNDERDILATED.
7.	Discard the used mounting and delivery system as biomedical waste (Bio-Hazard), incinerate the same and then handover to the biomedical waste handling department.

## 10. Disclaimer of Warranty and Limitation of Remedy

There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Translumina GmbH product(s) described in this publication. Under no circumstances Translumina GmbH shall be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law. No person has the authority to bind Translumina GmbH to any representation or warranty except as specifically set forth herein.

Descriptions or specifications in Translumina GmbH printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties. Translumina GmbH will not be responsible for any direct, incidental, or consequential damages resulting from reuse of the product.

Manufactured by:

**translumina GmbH** 

Neue Rottenburger Strasse 50

D-72379 Hechingen

Phone: ++ 49 74 71 98 94 - 0

Fax: ++ 49 74 71 98 94 - 3 80

E-mail: info@translumina.de

Internet: www.translumina.de

© 2020 translumina GmbH

All rights reserved



© YUKON® and translumina® are registered trademarks of translumina GmbH

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© YUKON® e translumina® são marcas comerciais registadas da translumina GmbH