


| | | | | |
|--|---|-------------------------------------|---|--|
| Product/Family: Yukon® Chrome PC Sirolimus Eluting Coronary Stent System | If applicable, subsystem / component: N/A | Manufacturer Translumina GmbH | | |
| Development project: Drug eluting CoCr-Stent | Project number: E060 | Delimitation, if applicable: N/A | | |
| Version number: 002 | Document purpose / Classification: TD according to MDR | Released on 2024-06-20 | Released by (area and name) Dr. Michael Stöver | Signature  |

1 Device Description and Specification (including Variants and Accessories)

1.1 Device Description and Specification

1.1.1) Technical specification

This chapter lists all the technical specifications of the Yukon® Chrome PC Sirolimus Eluting Coronary Stent System product. The technical specifications include the details of the dimensions, performance attributes, coating, specifications that relate to accessories and further characteristics of the sterile product.

Dimensions/ Design of the Yukon® Chrome PC Sirolimus Eluting Coronary Stent System product

This subsection lists the dimensions and design of the stent delivery system and stent of the Yukon® Chrome PC Sirolimus Eluting Coronary Stent System product.

The table below lists the specifications of the stent delivery system (SDS). The SDS has the same design for all sizes / variants of the Yukon® Chrome PC Sirolimus Eluting Coronary Stent System and differs only in the balloon length / balloon diameter. This difference also affects some other specified dimensions, therefore the range from smallest size to largest size is given for some attributes in the table below.

In addition, the table indicates whether the specifications are made available to the user (e.g. in brochures).

| Attribute | Specification | Available to user [YES / NO] |
|--|--|---------------------------------|
| Catheter usable length | 1360 mm to 1480 mm | YES |
| Marker for guiding catheter | | NO |
| a) from tip to distal marker | 909 mm ... 947 mm | |
| b) from tip to proximal marker | 1009 mm ... 1047 mm | |
| Length of guide wire lumen | 257 mm ... 295 mm | NO |
| X-ray marker distance | 8,5 mm ... 43 mm (for individual sizes s. product drawings) | NO |
| Crossing profile | 0.94 ... 0.96 mm (Ø 2.0 mm ... 2.5 mm; SV) 0.97 ... 1.20 mm (Ø 2.75 mm ... 4.0 mm; MV) Tolerance: + 5% | YES |
| Balloon working length | 9 mm ... 42 mm | NO |
| Proximal shaft diameter (Hypotube, PTFE coated) | 0.65 mm (1.9F) | YES |
| Distal shaft diameter | 0.87 mm to 0,91 mm (2.7F) | YES |
| Tip length | 3.5 mm ... 4.5 mm | NO |
| Tip entry profile | ≤ 0,49 mm | NO |
| Tip design | soft tip | YES |
| Balloon folding | trifold (memory) | NO |
| Hub | female luer lock catheter connector | NO |

Two stent designs of cobalt-chromium alloy L605 are used for the Yukon® Chrome PC Sirolimus Eluting Coronary Stent System:

- small vessel (SV) stent for diameter Ø 2.0 mm to Ø 2.5 mm
- medium vessel (MV) stent for diameter Ø 2.75 mm to Ø 4.0 mm

Within Small Vessel product range (2.0 mm to 2.5 mm) and within Medium Vessel product range (2.75 mm to 4.0 mm) there is no difference of stent design.

The table below lists the specifications of the two stent designs. The stent designs (SV & MV) are characterized by a homogeneous pattern of rings and links. Therefore, the only difference between shorter and longer stents is the number of identical rings. This difference in the number of identical rings also affects some other specified dimensions, therefore the range from smallest size to largest size is given for some attributes in the table below.

In addition, the table indicates whether the specifications are made available to the user (e.g. in brochures).

| Attribute | Specification | Available to user [YES / NO] |
|---|--|------------------------------|
| Strut thickness | 68 µm (SV) 79 µm (MV) | YES |
| Strut width (straight section) | 78 µm (SV) 88 µm (MV) | NO |
| Strut width (bow) | 88 µm (SV) 108 µm (MV) | NO |
| Strut width (connector) | 75 µm (SV) 76 µm (MV) | NO |
| Relative Stent surface (crimped) | 32% (Ø 2.0 mm) ... 33% (Ø 4.0 mm) | NO |
| Relative Stent surface (expanded) | 11.9 ... 14.9% (SV) 9.1 ... 13.3% (MV) | YES |
| Stent weight | 4.9 mg ... 19.0 mg (SV) 6.9 mg ... 32.9 mg (MV) | NO |
| Smallest approx. mesh diameter (for side branch access) | Ø 0.96 mm (SV) Ø 1.11 mm (MV) | NO |
| Cell circumference | 13.3 mm (SV) 18.4 mm (MV) | NO |
| Roughness | Homogenous microroughness by microporous surface ($R_{zmax} = 10 \mu m$) | NO |

Performance attributes of the Yukon® Chrome PC Sirolimus Eluting Coronary Stent System product

The table below lists the mechanical performance attributes of the Yukon® Chrome PC Sirolimus Eluting Coronary Stent System product. The specifications cover all test criteria applicable to the Yukon® Chrome PC Sirolimus Eluting Coronary Stent System by DIN EN ISO 25539-2.

In addition, the table indicates whether the specifications are made available to the user (e.g. in brochures).

| Attribute | Specification | Available to user [YES / NO] |
|---|--|------------------------------|
| Balloon deflation time | ≤ 30 sec. | YES |
| Balloon nominal pressure | 11 bar | YES |
| Balloon rated burst pressure | 16 bar | YES |
| Balloon fatigue / Multiple Inflation | 10 cycles at RBP without leak | NO |
| Dogboning | ≤ 5% | NO |
| Stent dislodgement force | ≥ 1N | NO |
| Tensile bond strength | Acc. To 10555-4 | NO |
| Torsional bond strength | Pass deflation spec. at 5 rotations | NO |
| Visibility X-ray | Pass acc. to 25539-2 | NO |
| Magnetic resonance imaging safety (MRI) | Pass acc. to 25539-2 | YES |
| Fatigue and durability of stent | Pass acc. to 25539-2 | NO |
| Crush resistance of stent (radial applied load) | ≥ 3 N | NO |
| Kink resistance (flexibility) of stent | Pass with Radius 5 mm (SV) / 10 mm (MV) acc. to 25539-2 | NO |
| Stent diameter to balloon inflation pressure (Compliance) | According to Compliance Chart (with tolerance +10% / -10%) | YES |
| Stent length (Foreshortening) | ≤ 5% | NO |
| Recoil | ≤ 7% | NO |

Specification of Coating of the Yukon® Chrome PC Sirolimus Eluting Coronary Stent System product

The coating solution is a mixture of:

- Poly(D,L-lactide) (Resomer R 202 S)
- Rapamycin (Sirolimus)
- Ethylacetate

with a ratio of 1:1, which is equivalent to a solution concentration of 4 mg Rapamycin and 4 mg Polylactide per 1 ml Ethylacetate.

As described in a separate process description coating is not applied to stent only, but to the stent, which already is crimped onto the balloon. The nominal Sirolimus content as specified in the product description of Yukon® Chrome PC therefore represents the overall amount of Rapamycin applied to stent and balloon.

The average ratio of the amount of Rapamycin on the stent is 48 % (T-CMG3016PC (Yukon® Chrome PC Sirolimus Eluting Coronary Stent System, \varnothing_{nom} : 3.0 mm, L_{nom} : 16 mm)) of the total amount of Rapamycin applied to stent and balloon. According to the specification the nominal Sirolimus content is 200 µg (T-CMG3016PC).

The stent is coated only on the abluminal surface.

The nominal amount of Rapamycin on stent surface is:

- 1.9 µg / mm² (stent surface)
- 6.2 µg / mm (stent length)

Specifications related to accessories

The table below lists the specifications related to accessories of the Yukon® Chrome PC Sirolimus Eluting Coronary Stent System product.

In addition, the table indicates whether the specifications are made available to the user (e.g. in brochures).

| Description | Specification | Available to user [YES / NO] |
|------------------------------|------------------------------------|------------------------------|
| Recommended guide wire | 0.36 mm (0.014") | YES |
| Recommended guiding catheter | 5 F | YES |
| Luer connection | Acc. To standard ISO 80369-1 to -3 | NO |

Further characteristics of the Yukon® Chrome PC Sirolimus Eluting Coronary Stent System product

The table below lists further characteristics specified for the Yukon® Chrome PC Sirolimus Eluting Coronary Stent System product.

In addition, the table indicates whether the specifications are made available to the user (e.g. in brochures).

| Description | Specification | Available to user [YES / NO] |
|-----------------------------|--|------------------------------|
| Number of applications | Single-use | YES |
| Particle burden | Specification on finished product incl. coating particles: >10 µm : max. 12000 >25 µm : max. 1200 >50 µm : max. 120 Specification for product without the coating: >10 µm : max. 6000 >25 µm : max. 600 >50 µm : max. 60 removable particles visible with the bare eye are not allowed | NO |
| Sterility | Sterile product | YES |
| Sterilization type | EO (ethylene oxide gas sterilization) | YES |
| Product shelf-life | 24 months | YES |
| Storage conditions | 8 – 25 °C (short term 40°C) | YES |
| Application time (stent) | > 30 days (permanent implant) | NO |
| Application time (catheter) | < 24 h | NO |

Revision History

| Rev.No. | Date | Created by | Description |
|---------|------------|------------|--|
| 001 | 2023-10-12 | PHA | New document for TD according to MDR |
| 002 | 2024-06-20 | MST | Adaption for inclusion of Natec variants |