EN	Instructions for use	DEB PTCA Balloon Dilatation Catheter
FR	Mode d'emploi	DEB Cathéter de dilatation pour ACTP
DE	Gebrauchsanweisung	DEB PTCA-Balldilatationskatheter (PTCA)
IT	Istruzioni per l'uso	DEB Catetere dilatatore à PTCA
NL	Gebruiksaanwijzing	DEB PTCA Ballon dilatationskatheter
ES	Instrucciones de Uso	DEB Catéter de dilatación (ACTP)
CS	Návod k použití	Balónkový dilatační katetr (PTCA)
EL	Οδηγίες χρήσης	DEB μπαλόνι διαστολής (PTCA)
PL	Instrukcja stosowania	DEB Balonowy cewnik do PTCA
PT	Instruções de Uso	Cateter Balão de Dilatação (PTCA)
TR	Kullanım talimatlarına	DEB PTCA Balon Dilatasyon Kateteri

**PROTÉGÉ**

DEB CATHETER

**PROTÉGÉ NC**

NON-COMPLIANT DEB

**BLUEMEDICAL**

*A Translumina Group Company*

Drug Eluting Percutaneous Transluminal Coronary Angioplasty Balloon Dilatation Catheter.

### Description

The distal section is dual lumen, one lumen permits the inflation of the balloon, the second lumen permits the use of a 0.014" (0.36 mm) guide wire, to facilitate the advance of the catheter to and through the stenosis. The proximal section of the catheter is a single lumen, stainless steel hypotube with a single luer port for connection to the inflation/deflation device. The balloon includes one or two radiopaque markers to aid in the positioning of the balloon within the stenosis. Marks located on the proximal portion of the catheter shaft indicates the exit of the balloon catheter tip out of the guide catheter (one at 90cm and one at 100cm). An active component, Paclitaxel, is incorporated in a coating on the balloon surface. On inflation of the balloon the coated surface comes in contact with the vessel and the incorporated Paclitaxel is delivered to the vessel wall. The total load of Paclitaxel is dependant on diameter and length of the Protégé balloon. The coated part of the balloon holds 3µg Paclitaxel per mm<sup>2</sup>.

*Drug load*<sub>x</sub> = π \* ϕ<sub>balloon</sub> \* length<sub>balloon</sub> \* 3

*Drug release*<sub>max</sub> = π \* ϕ<sub>balloon</sub> \* length<sub>balloon</sub> \* 0,27

Balloon size	Drug load	Max drug release
2.00x10mm	188 µG	17 µG
4.00x30mm	1131 µG	102 µG

### Intended Purpose / Intended Use

The DEB PTCA catheter is intended for balloon dilatation of the stenotic portion of a coronary artery or bypass graft stenosis for the purpose of improving myocardial reperfusion in patients during a percutaneous coronary interventional (PCI) performed in a blood vessel. The DEB PTCA catheter is intended for transient use (≤60 minutes).

### Indications

- Indicated for balloon dilatation of the stenotic portion of a coronary artery or bypass graft stenosis, including in–stent restenosis (ISR), for the purpose of improving myocardial reperfusion in de novo lesions and small vessels (SVD)
- Pre– and post balloon dilatation during coronary stent procedures

### Contra–Indications

The use of the balloon dilatation catheter is contraindicated for:

- Intolerance or allergy to Paclitaxel
- Unprotected left main coronary artery
- Coronary artery spasm in the absence of a significant stenosis
- Patients that are (potentially) pregnant

### Patient Target Groups

The DEB PTCA catheter is intended to be used in patients suffering from stenotic lesions in the coronary vascular system who require a PCI.

### Intended User(s)

The target users are medical professionals (e.g. Interventional Cardiologists, cardiovascular technicians) who perform and assist in catheterisation procedures in a clinical laboratory setting, during percutaneous intervention (e.g. PTCA, coronary stent placement).

Use of the product is restricted to healthcare professionals in an interventional suite with angiography equipment, with personnel with relevant and adequate training and who are familiar with the possible/conceivable complications.

### Clinical Benefits

Balloon dilatation of the target lesion with concomitant release to the vessel wall of the antiproliferative drug, paclitaxel, with the benefit of a lower risk of target lesion restenosis (TLR) in the mid– (>6 months) and longer term (> 12 months) when compared to the TLR outcomes of non–drug coated balloon dilatation (conventional PTCA/POBA) as reported in the literature. This means for the patient a lower risk of a reintervention due to TLR.

### Warnings

- For single patient, single procedure use only. Do NOT resterilize and/or reuse as this can potentially result in compromised device performance and increased risk of inappropriate resterilization and cross contamination.
- To reduce the potential vessel damage, the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis.
- PTCA in patients who are acceptable candidates for coronary artery bypass graft surgery requires careful consideration. This includes possible hemodynamic support during PTCA as treatment of this patient population carries special risk.
- When the balloon dilatation catheter is exposed to the vascular system, it should be manipulated while under high quality fluoroscopic guidance. Do not advance or retract the catheter unless the balloon is fully deflated under vacuum. If resistance is met during manipulation, determine the cause of the resistance before proceeding.
- Balloon Pressure: do not exceed the rated burst pressure (RBP) indicated on the packaging. The rated burst pressure (RBP) is based on results of in vitro testing. At least 99.9% of the balloons (with a 95% confidence) will not burst at or below their rated burst pressure (RBP).
- Use of a pressure–monitoring device is recommended to prevent over pressurization.
- PTCA should only be performed at hospitals where emergency coronary bypass surgery can be quickly performed in the event of a potentially injurious or life–threatening complication.
- Use only the recommended balloon inflation medium. Never use air or any gaseous medium to inflate the balloon.
- Use the catheter prior to “Use Before” (“expiry”) date specified on the package.
- Do not use if packaging is damaged or unintentionally opened prior to use.
- Treatment of vessels with moderate or heavy lesion calcification is associated with decreased success rates up to 60–85% and increases the risk of acute closure, vessel trauma, balloon entrapment and associated complications (<1%).
- Paclitaxel is in potential a genotoxic (in particular clastogenic) agent based upon its pharmacodynamic mechanism of action, which is interference with the microtubule disassembly. The relevance of this specific mechanism of genotoxicity for human carcinogenity risk is currently not known.
- Breastfeeding should be withheld for at least 6 days after treatment.

### Interactions with other drugs

The total load of Paclitaxel is only a fraction of less than 0.03% of the volume generally applied in chemotherapies. The risk of interactions with other active agents is therefore high improbable. Nonetheless caution should be exercised when concurrently administering known CYP3A4 and/or CYP2C8 substrates (like cyclosporine, lovastatin, midazolam, ondansetron, terfenadine) or drugs with high PPB (like sulfonureas, coumarin type anti–coagulants, digitoxin, salicylic acid, sulfonamides) the specific instructions for use for these active agents should be consulted in addition.

### Precautions

- Prior to angioplasty procedure, the catheter should be examined to verify functionality and ensure that its size and shape are suitable for the specific procedure for which it will be used.
- Only physicians trained in the performance of percutaneous transluminal coronary angioplasty should use the catheter.
- Direct touching of the balloon, wiping of balloon surface or contact with organic solvents, e.g. alcohol, should be strictly avoided, since this can cause delamination of the balloon coating.
- During the procedure, appropriate anticoagulant and coronary vasodilator therapy must be provided to the patient as needed. Anticoagulant therapy should be continued for a period of time to be determined by the physician after the procedure.

### Potential complications and adverse effects

The following complications, among others, can occur by performing of a Percutaneous Transluminal Coronary Angioplasty

- Death.
- Acute myocardial infarction.
- Total occlusion of the coronary artery or bypass graft.
- Coronary vessel dissection, perforation, rupture or injury.
- Restenosis of the dilated vessel.
- Haemorrhage or haematoma.
- Unstable angina.
- Arrhythmias, including ventricular fibrillation.
- Drug reactions, allergic reaction to contrast medium.
- Hypo/hypertension.
- Infection.
- Coronary artery spasm.
- Arteriovenous fistula.
- Embolism.

### Pharmaceutical information on included drug

Minor traces of Paclitaxel in the blood plasma appear less relevant than for a systemic treatment. However, side effects cannot be completely ruled out.

- Abnormal liver enzyme values.
- Allergic or immunological reaction to the drug or similar agents.
- Alopecia.
- Anaemia.
- Disorders of the heart conduction system.
- Gastro–intestinal tract impairment.
- Haematological dyscrasia (incl. leukopenia, neutropenia, thrombocytopenia).
- Histological changes in the vascular wall, incl. inflammation, cell damage or necrosis.
- Myalgia / arthralgia.
- Peripheral neuropathy.
- Pseudomembranous colitis.

### Inspection prior to use

Prior to use of the balloon dilatation catheter in PTCA, all equipment to be used for the procedure, including the dilatation catheter should be carefully examined for defects. Examine the dilatation catheter for bends, kinks or other damage. Do not use any defective equipment. Prepare equipment to be used following manufacturer’s instruction or standard procedure.

### Preparation of the Dilatation catheter

Complete the following steps to prepare the catheter:

- Remove dilatation catheter from the package.
- Remove catheter from hoop.
- Remove protective sheath with the stylet from balloon catheter.
- Prepare an inflation device with the recommended contrast medium according manufacturer’s instruction.
- To evacuate air from the balloon segment, the following procedures should be followed.
  - Fill a 10cc syringe with approximately 2cc 50/50–mixture contrast/saline.
  - Connect the 10cc syringe to the balloon dilatation catheter, hold the nozzle downwards and pull the plunger to create a negative pressure (hold this for at least 15 seconds).
  - Slowly release plunger to neutral (with nozzle of syringe pointed downwards), allowing contrast to fill the shaft of the balloon dilatation catheter.
  - Disconnect syringe from balloon dilatation catheter and remove air out of syringe.
  - Repeat the above points until no more bubbles appear in the syringe.
  - Disconnect syringe from the balloon catheter and connect the inflation device, with a fluid–fluid connection, to the balloon catheter.
  - Pull the plunger of the inflation device to create a negative pressure.

### Warnings

- An incorrectly prepared balloon catheter may prolong the deflation time.
- Loss of vacuum or a continuous stream of air entering the syringe upon aspiration indicates the presence of a leak in the system (inspect the system for integrity).
- Do not torque the balloon catheter shaft.

### Insertion technique

- First advance the guiding catheter and connect to Y–connector (minimal 0.056" / 1.42 mm ID) and cross the stenosis with a 0.014"/ 0.36 mm (or smaller) guide wire according to the manufacturer's instructions.
- Wipe the extending guide wire with a sponge or gauze to remove residuals of blood or contrast.
- Load the guide wire into the distal tip of the balloon dilatation catheter guide wire lumen and ensure that the guide wire exits at the transition point 25 cm proximal of the distal tip.
- Advance the balloon dilatation catheter over the guide wire until it approaches the Y–connector.
- Loosen the valve of the Y–connector and insert the balloon dilatation catheter for minimal 30 cm while maintaining the guide wire position.
- Close the Y–connector to create a seal around the balloon dilatation catheter but ensure that movement is not restricted.
- Advance the balloon to the stenosis under fluoroscopy.
- Inflate the balloon while monitoring the inflation pressure, the balloon position, the ECG, blood pressure and other vital conditions of the patient.
- After dilatation of the stenosis deflate the balloon, check deflation under fluoroscopy.
- Pull the balloon dilatation catheter back, into guiding catheter, while maintaining the guide wire position.
- Perform a control angiography.
- Continue the procedure according to the accepted PTCA protocol.

### Catheter exchange procedure

- Loosen the haemostatic valve.
- Hold the guide wire and haemostatic valve in one hand, while holding the balloon shaft in the other hand.
- Maintain guide wire position in the coronary artery by holding the wire stationary and pull the balloon dilatation catheter out of the guiding catheter while monitoring wire position under fluoroscopy.
- Withdraw the deflated balloon dilatation catheter until the transition point in the guide wire lumen is reached. Carefully pull the flexible, distal portion of the balloon catheter out of the haemostatic valve while maintaining the guide wire position across the lesion.
- Prepare the next balloon catheter to be used, as previously described.
- Backload the new balloon catheter onto the guide wire as previously described and continue the procedure accordingly.

### Warning

After use, the equipment used in angioplasty, may be a potential biohazard. Handle carefully and dispose of in accordance with hospital policy and appropriate, applicable federal regulations.

### Notice to the user

In case of any serious incident that has occurred in relation to the device, this should be reported to the manufacturer and the competent authority of the Member State in which the user and/ or patient is established.

**Link to the Summary of Safety and Clinical Performance (SSCP)**
The Summary of the Safety and Clinical Performance (SSCP) is available in the European database on medical devices (Eudamed).
The DEB PTCA catheter has Basic UDI-DI **872063420BMD–DEB6S**.
**https://ec.europa.eu/tools/eudamed**

### Warranty

Blue Medical Devices B.V., (BLUE), warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excluded all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties or merchantability of fitness. Handling storage, cleaning and sterilization of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond BLUE control directly affect the device and the results obtained from its use. BLUE's obligation under this warranty is limited to the repair or replacement of this device and BLUE shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this device. BLUE neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. BLUE assumes no liability with respect to devices reused, reprocessed or re–sterilized and makes no warranties expressed or implied, included but not limited to merchantability or fitness for intended use, with respect to such device.

**Manufactured by:**
Blue Medical Devices B.V.
Panovenweg 7
5708 HR Helmond
The Netherlands
Phone: +31 (0) 492 588900

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Medikamentenfreisetzender perkutaner transluminaler Koronarangioplastie-Ballon-Dilatationskatheter.

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© Blue Medical Devices B.V.
Blue Medical Devices B.V. is an