

Trans-ID

Balloon Inflation Device - Gun Type



Maximum Pressure of 30 atm

20 ml clear syringe barrel

Trans-ID self locking mechanism for easy locking

Precision manometer

Inflation- Easy to inflate even at Higher pressure

Deflation- Immediate quick trigger pressure release



(With Kit)



Sr. No	The Name Of The Medical Product	Additional features / Requirements / Technical Specification		Expiration Date
01	Inflation Device	Maximum Pressure:	30 ATM	Expiration date medical device on the delivery date should be at least 75% (seventy-five) from general term storage or at least 15 (fifteen) months
	Inflatable Device	Full Kit:	Pump, 3-way crane, Y-connector, input device conductor, rotating device	
		Display:	Located at an angle of 30 °, with a fluorescent background. The presence of modifications from analog and digital LED displays	
		Flexible Tube With 3-way Valve:	Length not less 33 cm	

Sr.no	PRODUCT DESCRIPTION	PRODUCT CODE
01	Trans-ID Inflation Device without kit	TRID30G
02	Trans-ID Inflation Device with push Pull kit	TRID30G1
03	Trans-ID Inflation Device with Screw Type Kit	TRID30G2

Trans-ID

Balloon Inflation Device


Maximum Pressure of 30 atm

20 ml clear syringe barrel
with smooth moving plunger

Trans-ID self locking mechanism
for easy locking as well as gradual
and rapid pressure

Ergonomic handle

Precision manometer



Reference No.	Pressure (atm)	PRODUCT CODE
TRID30	30	20

Trans-ID

Balloon Inflation Device - Gun Type

Maximum Pressure of 30 atm

20 ml clear syringe barrel

Trans-ID self locking mechanism for easy locking

Precision manometer

Inflation- Easy to inflate even at Higher pressure

Deflation- Immediate quick trigger pressure release




Sr.no	PRODUCT DESCRIPTION	PRODUCT CODE
01	Trans-ID Inflation Device without kit	TRID30G
02	Trans-ID Inflation Device with push Pull kit	TRID30G1
03	Trans-ID Inflation Device with Screw Type Kit	TRID30G2

Trans-ID

Balloon Inflation Device - Normal Type

Accurate Pressure Control
50% Faster For Inflation

Easy for Deflation
Latex-Free



Sr. No	The Name Of The Medical Product	Additional features / Requirements / Technical Specification		Expiration Date
01	Inflation Device	Maximum Pressure:	30 ATM	Expiration date medical device on the delivery date should be at least 75% (seventy-five) from general term storage or at least 15 (fifteen) months
	Inflatable Device	Full Kit:	Pump, 3-way crane, Y-connector, input device conductor, rotating device	
		Display:	Located at an angle of 30 °, with a fluorescent background. The presence of modifications from analog and digital LED displays	
		Flexible Tube With 3-way Valve:	Length not less 33 cm	

Sr.no	PRODUCT DESCRIPTION	PRODUCT CODE
01	Trans-ID (Normal)	TRID30A