




trans-ASPIRATOR

Aspiration Catheter Kit



LARGE ASPIRATION LUMEN

Constant, high performance Aspiration offers Quick & complete Aspiration

OPTIMIZED TIP DESIGN WITH RADIO OPAQUE PLATINUM IRIIDIUM MARKER BAND

Smooth Tip design enhancing deliverability

Radio Opaque Platinum Iridium Marker Band ensures reliable fluoroscopic visibility

BRAIDED SHAFT WITH HYDROPHILIC COATING

Braided Shaft ensures excellent pushability and offers kink resistance

30 CM hydrophilic coating from the tip ensures optimum flexibility and Trackability

WITH STYLET (REMOVABLE CORE WIRE)

Removable core wire enhances Pushability and Anti kinking feature.

Catheter length: 145 cm

Rapid exchange length: 25 cm

Summary of Technical Specifications

Usable Length	145 cm
Guiding Catheter Compatibility	For 6F Min. ID = 0.071" For 7F Min. ID = 0.082"
Distal Tip Hole Length	3mm for 6F
Rapid Exchange Segment	25 cm
Guide wire compatibility	0.014"
Coating	Hydrophilic